

Treadmill Gait Lab (Web App)

Step & Stride Metrics (Side View, One Runner)

What you'll measure (treadmill running)

- Step time (ms): time between alternating foot strikes ($R \rightarrow L$ or $L \rightarrow R$).
- Stride time (ms): time between same-foot strikes ($R \rightarrow R$ or $L \rightarrow L$).
- Stride frequency (Hz): $1 / \text{stride time (sec)} = 1 / (\text{stride time (ms)} / 1000)$.
- Step length (m) (treadmill-derived): $\text{speed (m/s)} \times \text{step time (sec)}$.
- Stride length (m) (treadmill-derived): $\text{speed (m/s)} \times \text{stride time (sec)}$.

Note: On a treadmill, the runner stays in place, so lengths are derived from belt speed \times time (not from forward displacement on the floor).

Equipment & setup

- Phone/laptop with camera + browser (Chrome/Safari recommended).
- Stable support (tripod / propped against a stable object).
- Known treadmill speed (km/h or m/s).

Camera position (critical)

- Side view only (sagittal plane).
- Keep hips-to-feet visible at all times.
- Avoid treadmill rails blocking ankles/heels.
- Good lighting; avoid strong backlight (bright window behind runner).

Procedure (Live Camera)

1. Open the web app link and allow Camera access.
2. Enter treadmill speed (km/h or m/s).
3. Start treadmill at a steady pace; runner moves naturally.
4. Tap Start Live Camera.
5. Wait until the table records 10 steps (auto-stops), or press Stop.
6. Tap Download CSV to save results.

Procedure (Upload Recorded Video)

7. Record a side-view treadmill clip (10–15 seconds is enough).
8. Open the app \rightarrow Upload video \rightarrow click Process Uploaded Video.
9. Let it run until 10 steps are captured \rightarrow download CSV.

Quality checks (before you trust your data)

- Tracking quality should show Good or Medium.
- Feet should be visible and not blurred badly.
- If steps are missed or double-counted: move camera slightly back, raise camera height slightly, improve lighting, or reduce occlusion from rails.
- Facing indicator: if it stays Unknown, adjust camera angle/lighting so one body side is clearly visible.